

SPARNET – Spartan Data Network for Real-Time Physiological Status Monitoring

Reed W. Hoyt, Ph.D.

Biophysics and Biomedical Modeling Division

US Army Research Institute of Environmental Medicine (USARIEM)

Natick, MA 01760-5007

email: reed.hoyt@us.army.mil

U.S. Army Telemedicine Partnership Series 2008: “Personal Health Monitoring”

Seattle, Washington

April 5, 2008



Report Documentation Page				Form Approved OMB No. 0704-0188	
Public reporting burden for the collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to a penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number.					
1. REPORT DATE 05 APR 2008		2. REPORT TYPE		3. DATES COVERED 00-00-2008 to 00-00-2008	
4. TITLE AND SUBTITLE SPARNET - Spartan Data Network for Real-Time Physiological Status Monitoring				5a. CONTRACT NUMBER	
				5b. GRANT NUMBER	
				5c. PROGRAM ELEMENT NUMBER	
6. AUTHOR(S)				5d. PROJECT NUMBER	
				5e. TASK NUMBER	
				5f. WORK UNIT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) US Army Research Institute of Environmental Medicine (USARIEM), Biophysics and Biomedical Modeling Division, Natick, MA, 01760-5007				8. PERFORMING ORGANIZATION REPORT NUMBER	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)				10. SPONSOR/MONITOR'S ACRONYM(S)	
				11. SPONSOR/MONITOR'S REPORT NUMBER(S)	
12. DISTRIBUTION/AVAILABILITY STATEMENT Approved for public release; distribution unlimited					
13. SUPPLEMENTARY NOTES Army Telemedicine Partnerships Series 2008: ?Personal Health Monitoring?, Seattle, WA, 5 Apr 2008					
14. ABSTRACT					
15. SUBJECT TERMS					
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT Same as Report (SAR)	18. NUMBER OF PAGES 27	19a. NAME OF RESPONSIBLE PERSON
a. REPORT unclassified	b. ABSTRACT unclassified	c. THIS PAGE unclassified			

Outline

1. Overview of USARIEM
2. Need for real-time physiological monitoring
 - Ranger Training
 - CST-WMD
3. Overview of Spartan data networks (SPARNET)
 - Inductive-Personal Area Network
 - Software-defined Squad Area Network
4. Conclusions



U.S. Army Research Institute Environmental Medicine (USARIEM)



- A subordinate laboratory of the U.S. Army Medical Research and Materiel Command
- Mission: conduct basic and applied research to determine how extreme heat, severe cold, high terrestrial altitude, occupational tasks, physical training, deployment operations, and nutritional factors affect the health and performance of military personnel



Physiological Monitoring and Predictive Modeling

Develop wearable physiological sensor systems that collect, organize and interpret data on Soldier health status



Develop thermoregulatory models and decision aids to predict risk and performance decrements for Soldiers under the stresses of heat, cold and high altitude


Real-time physiological monitoring applications



Ranger Training Brigade, Fort Benning, GA

Land Navigation Training Area

Mitigate risk of injury and separation by tracking Ranger student physiological status and physical location during higher risk training activities.



5 Mile Run - 8 min/mile



ypoints
Routes
GPS
Tracks(2)Nan
ACT
ACT

Core Temp (°F)

104
103
102
101
100
99
98
97
96

40 min Run

5 Mile Run (about 0545-0700)

Test volunteers:

Age: 23 ± 3 yr; Wt: 167 ± 30 lbs, Body fat: 16 ± 3 %, Est. VO_2 max: 55 ± 3 ml/kg/min

Weather:

Tair = $71 \pm 1^\circ\text{F}$, Humidity = 95% RH


Civil Support Team–Weapons of Mass Destruction (CST-WMD)

MEDICAL MONITORING TELEMETRY SYSTEM


Proactive medical monitoring can decrease serious injuries and loss of operational strength/readiness by tracking changes and trends in vital signs/core temperatures and taking action before they become symptomatic and result in a medical emergency "man-down" scenario



Medical Monitoring Telemetry



Spartan data networks (SPARNET)



SPARNET Goal

Develop a state-of-the-art wireless squad and personal area network that enable existing physio-sensors, algorithms, and interfaces to form real-time medical and situational awareness products for foot Soldiers

- Use network monitoring of physiology and geolocation to:
 - Minimize casualties and improve casualty management
 - Improve mission planning and risk management
 - Increase mission training opportunities



SPARNET Guiding Principles

- Minimalism
 - Focus requirements to minimize risk/cost
- Flexibility and adaptability technologies
 - Open, easily-modified, “white box”
 - Facilitates response to specialized needs
- Focus on dismounted soldier needs
 - Establish close partnerships



Multidisciplinary “triple helix” Community Infrastructure Project

Government (Army)

- **Military customers, program leadership, science and engineering expertise** (textile prototypes, physiological monitoring and biomedical modeling, data base)
- **Regulatory and safety oversight** (RF safety, Army Spectrum Manager, IRB)
- **Funding:** SBIR/STTR, TATRC, MRMC

Small business


- **Engineering leadership and expertise:** Elintrix.com
 - Software defined radios (RF and inductive networks)
 - Situational Awareness application

Academia

- **Wearable antenna design:** Ohio State Univ. ElectroScience Laboratory
- **Power efficient radio and network operation:** Univ. of California San Diego
- **Metabolic rate algorithms:** Rice University




Personal Area Sensor Network



VitalSense® Integrated Physiological Monitoring System (MiniMitter.com)

- Ingestible temperature pill
- Skin temperature patch
- Push PAN (unidirectional)
- Near field RF
- Resistant to interference
- Data logger
- FDA 510k



Warfighter Physiological Status Monitoring (WPSM)



Vital Sign Detection System (Hidalgo, Ltd, UK)

Sensor node: ECG, respiration, actigraphy,
body orientation, skin temp



← **Core Temp** (MiniMitter)
Ingestible Thermometer Pill




← **Fluid Intake Monitor**



← **Sleep Watch**

Future → Inductive-PAN

- Free space magnetic induction
 - Inherently short range/near-field
 - Body Bubble
 - Spatial reuse of spectrum
 - Tissue transparent
 - Issue: No adequate COTS transceiver available
 - BioNet – ground-up nodes being developed (STTR)



SPARNET Squad Area Network (SAN)

LAND NAV

TRAINING/FTX

- LOST
- SEPARATED
- PANIC
- INJURY
- STATIONARY

Alert Mechanism

**Tactical Operations
Center (TOC)**

**Laptop
Computer**

**Backbone
Lightweight
Repeaters
(Data Only)**

Medics

GPS


PDA

Squad of Ranger Students

Ranger Instructor





How a "Cognizant Radio" Helps





SAN Radio and Antenna

- Software-Defined-Radio (SRD)
- Near-Field-Communication Hub (forms links to sensors)
- Enclosure, seal, circuit card assemblies, lower cover, batteries, battery cover
- Engineering unit LxWxH: 6"x3.7"x1"




SAN Body-Worn Antenna

Front Position: Double Null





Shoulder Position: Reduced Nulls



Risk-Reduction by Extensive Analysis, Modeling and Simulation



Angle Positioning



Best location: $\theta = 165^\circ$



Squad Area Network (SAN) Elements

Wearable
antenna

SAN module

AA
Batteries



SPARNET Initial Objectives

- **Battery life: >18 hours (TBD)**
- **Radio Unit (SDR with cognizant-radio capabilities)**
 - Flexible modulation options
 - Prototype size: 6" x 3.7" x 1" (~15 cm x 9 cm x 2.5 cm)
 - Weight: ~1kg
 - Inductive network links
 - Target single-hop range: >300 meters
 - Target soldier-911 range: ~1000 meters
 - Prototype operating band: 225 MHz – 380 MHz
 - Frequency de-confliction: user-selected keep-out bands
- **Antenna**
 - Body-worn, integrated into hydration unit
 - Link to external repeaters for increased performance



SPARNET IN SUMMARY...

- Risk-managed to address significant unmet needs
- Minimalist, adaptable, transparent, advanced technology
- Scalable and cost-effective
- Improves training effectiveness and safety
 - **Situational and medical awareness**
 - **Casualty avoidance**
 - **Casualty management**
- Enables innovative, cutting-edge physiologic research
- Enables cutting-edge network science research





Acknowledgement

- Mark Buller
- Drew Barnett
- Stephen Mullen
- Chi-Chih Chen
- Ranger Training Brigade
- TATRC



Citations of commercial organizations and trade names in this report do not constitute an official Department of the Army endorsement or approval of the products or services of these organizations.

Opinions, interpretations, conclusions, and recommendations contained herein are those of the author and are not necessarily endorsed by the U.S. Army

